ReFa

ReFa RHYTHM

GUIDE BOOK

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GUIDE BOOK

English

Thank you for purchasing ReFa RHYTHM.

ReFa represents a sincere desire to see women harness the beauty that flows from within and without them——.

It is a brand that expresses a wish to see you delight in discovering your true radiance with just a minimum of effort.

It also symbolizes an intention to produce a spontaneous smile on your face as you enjoy a sense of fulfillment from simply using our product.

ReFa is here for your everyday pleasure.

■ Pushing and rolling action

This device features a structure with 25 protruding projections on the left and right rollers and is designed after factoring in the size and intervals required to loosen the muscles of a small and thin face. The roller facilitates blood circulation by pressing the conical projections on your skin, producing a relaxing effect. You can also use the roller to care for your skin by gently kneading the desired areas ^{*1}.

■ Microcurrent ^{*2}

ReFa rollers emit a microcurrent by taking in light through a built-in solar-panel.

■ Platinum Coating (Roller)

■ Waterproof Construction

Conforms with JIS standards (equivalent to IPX7).


*1 The effects of finger-pressure therapy (improved blood circulation, relief of muscle pain, loosening of muscle tightness) cannot be achieved unless you press the conical projections on the rollers against your skin.

*2 Microcurrents have no effect on the improvement of blood circulation, relief of muscle pain, or loosening of muscle tightness.

Rolling Method

The strength in which the rollers pull up your skin can be adjusted by changing the angle at which you hold the handle.


<When Pushing the Rollers Forward>

Gradually raise the end of the handle to check the pull of the rollers.

*Lower the handle if you begin to feel pain.

<When Drawing the Rollers Back>

Gradually lower the end of the handle to check the pull of the rollers.

*Raise the handle if you begin to feel pain.

■ Applying the Rollers to your Skin

- (1) Grip the handle lightly and touch the rollers to your face.
*Avoid bending your wrist and hold the handle in a natural position.
- (2) Roll the rollers slowly by moving your arm up and down.

Refer to the arrows in the picture when rolling.

-  : Rolling back and forth
-  : Perform rolling by gradually increasing pressure from the narrow part of the arrows to the wide part.
-  : Lightly push the conical projections against your skin

How To Use

- Refer to the instructions in the following pages and use according to your preference.
- As a general rule, roll back and forth over each body part about 5 to 10 times.

Please be sure to read the instruction manual before use

- This product is intended for use on the **face only**. Never use this product on any parts of the body other than those directed in this instruction manual. **Doing so may result in internal bleeding and bruising, or other injury.**
- Start off by using the device gently for a short period of time while observing the condition of your skin. The conical projections on the rollers, if used with strong force or for a long time, **may cause internal bleeding and bruising, or other injury.**
- Depending on predisposition or your physical condition, internal bleeding, bruising, swelling, or itching may occur during or after use. Before using the device, be sure to test it on an inconspicuous part of the body to check that no abnormal reactions occur. If you experience any abnormal reactions, stop using the device and consult with your physician. *Some temporary redness may occur after usage. If this happens, use the device while paying attention to the condition of your skin.
- Persons with metal allergies, or who have conditions such as skin infections or dermatitis, should not use the product.
- Women in the early unstable stages of pregnancy or immediately after childbirth, or women who are pregnant and not in good health should not use this product.
- Remove any jewelry or accessories such as necklaces and piercings during use as they may get caught in the device.
- Tie your hair up or wear a hair band when using the device so your hair does not get caught.

Face

Loosen the stiffness in the face for a refreshing feel*.

Hold your chin in the rollers and roll them horizontally across the chin bit by bit.

Roll the device along the jawline towards the part below your ears by taking care that your chin is caught in the roller.

Roll along the cheekbones up to the sides of the ears. Push the projections of the rollers against the ○ portion.

Roll along the jaw bit by bit. At the end, you can use one side of the roller to knead the ○ portion.

* The effects of finger-pressure therapy (improved blood circulation, relief of muscle pain, loosening of muscle tightness) cannot be achieved unless you press the conical projections on the rollers against your skin.


Push the rollers against the portion above the eyebrows.

Close your eyes, stretch your forehead and roll bit-by-bit and in a zigzag manner. Finally, roll it across the hairline.

* The effects of finger-pressure therapy (improved blood circulation, relief of muscle pain, loosening of muscle tightness) cannot be achieved unless you press the conical projections on the rollers against your skin.

Main muscles

*The illustration uses generic names and does not imply any effect of the product.

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